

Laboratory Services

Important Test Announcement

FENTANYL, URINE

New Test Code(s): FENTB

Live Date: 10/13/2025

- **Methodology:** Enzyme immunoassay (EIA)
- **Performed:** Daily
- **Turnaround Time:** 1 Day
- **Specimen Required:** Urine collected in a sterile container, 10 mL (5 mL minimum)
 - **Specimen Preparation:** No special collection
 - **Storage/Transport/Temperature/Conditions:** Specimens can be transported to the lab at room temperature and are stable when stored at 4°C for up to 7 days after collection
 - **Unacceptable Conditions:** Samples other than urine, use of a urine preservative, collected in diaper or cotton ball
- **Comments:** **FENTB** will be added to all Urine Drugs of Abuse panels which include DOA, DOAT, DOAC, and DOATC, but can also be ordered as a stand-alone test.
 - All urine drug testing performed by this method is for screening purposes only. Unconfirmed results should only be used for medical purposes. False positive and erroneous results can occur due to cross reactivity.
 - **Reference Intervals:** Presumptive Negative

Please note, this is a screening test. If confirmatory testing is required, please order **XFENTB** on the earliest available urine sample. If necessary, consider ordering as an add-on to the screening sample.

If you have any additional questions about **Fentanyl, Urine**, please refer to the Laboratory Test Directory or call Client Services at 614-722-5477.