



## Laboratory Services

### Important Test Announcement

#### IBD PRECIS

**New Test Code: XPIDB**

**Live Date: 5/26/2026**

Effective May 26, 2026, IBD sgi Diagnostic (**XIBDB**) testing will be discontinued as it no longer offered by the reference laboratory. A new test, IBD PreciS (**XPIDB**) will be the replacement test. Sample collection requirements will change. Refer to the chart below and our Laboratory Test Directory for complete collection information.

	<b>New Test IBD PreciS XPIDB</b>	<b>Old Test IBD sgi Diagnostic XIBDB</b>
<b>Methodology</b>	High Performance Liquid Chromatography/ Tandem Mass Spectrometry (LC-MS/MS)	Multiple methods
<b>Performed</b>	Monday through Friday	Monday through Friday
<b>Turnaround Time</b>	4 days	4 days
<b>Specimen Required</b>	Gold Tube, (SST) with gel	3 mL Purple Tube, no EDTA  AND 6 mL Red Tube, no gel
<b>Specimen Volume</b>	3.5 mL	9 mL

<b>Specimen Preparation</b>	Centrifuge, remove serum and ship refrigerated.	Do not centrifuge purple EDTA tube  Transport to laboratory as soon as possible
<b>Storage/Transport/Temperature</b>	Refrigerated	Refrigerated
<b>Unacceptable Conditions</b>	None listed	Not all required tubes submitted
<b>Stability</b>	Whole Blood: Room temperature 3 days Whole Blood: Refrigerated 3 days  Serum: Room temperature 1 week Serum: Refrigerated 2 weeks	Whole blood: Red tube, (no anticoagulant), No Gel- Room temperature 3 days Whole blood: Red tube, (no anticoagulant), No Gel- Refrigerated 3 days  Whole blood: Purple tube (EDTA)-Room temperature 1 week Whole blood: Purple tube (EDTA)-Refrigerated 3 weeks  Serum: Room temperature 1 week Serum: Refrigerated 3 weeks

If you have any additional questions about **XPIBD**, please call Client Services at 614-722-5477.