



Please indicate the order of priority (1, 2 or 3) of the preparation type if you are requesting more than one preparation type for the same tissue specimen.

#### REQUEST INFORMATION: PREPARATION DETAILS

**FRESH PREP TYPE:**  Required  If available

Standard Fresh Preps:  RPMI  DMEM  Dry  PBS  Saline  RNALater  Other (please specify)  
 Slides-Touch Preps (# req'd )

**FROZEN PREP TYPE:**  Required  If available

Standard Freezing Methods:  Snap frozen (LN2 vapor phase)  OCT  
Slides:  H&E slides (# req'd )  unstained (# req'd )  Touch Prep (# req'd )  
Non-standard Freezing Methods:  Frz -20°C  Frz -80°C  Frz in isopentane  Dry Ice  
 Scroll/Ribbon (@ microns)  Other

**FIXED PREP TYPE:**  Required  If available

Standard formalin-fixed, paraffin embedded (FFPE)  
 Other (non-standard fixative, free floating, etc), specify:  
Slides:  H&E slides (# req'd )  unstained (# req'd )  
 Scroll/Ribbon (@ microns)  Other

#### SPECIMEN SIZE AND SHIPPING

Weight (preferred for solid tissue): Min Max  grams  mg

Dimensions: L W H

Volume (fluids): Min Max  mL  µL

Target # of Donors for this Request:

Saturday delivery:  Yes  No  If notified

Shipping Instructions:  Frozen  Refrigerated  Room temperature

Shipping Choices:  Ship Day of Procurement to Arrive Next Day (standard for fresh shipments)

Standard Overnight Shipment  Investigator pickup same day as procured