



Microbiology Swab Collection Chart

BACTERIAL/FUNGAL COLLECTIONS




COPAN Dual Culture Swab System with Liquid Stuart's Medium


CHLAMYDIA AND GONORRHOEAE COLLECTIONS BY AMPLIFIED DETECTION



APTIMA Unisex Swab Specimen Collection Kits for Endocervical and Male Urethral Swab Specimens




APTIMA Multitest Swab Specimen Collection Kit (vaginal swab specimens)




VIRAL COLLECTIONS



Rapid Flu – Puritan Sterile Foam-Tipped Applicator Swab



Rapid RSV – Puritan HydraFlock Applicator Swab



Viral Transport Media Collection Swab – Fisherbrand Polyester-Tipped Applicator Swab



Viral Transport Media Collection Swab – COPAN Nylon Flocked Mini-Tipped Flexible Plastic Applicator Swab



For more information regarding microbiology swab collections, please call
(800) 934-6575 or visit NationwideChildrens.org/Lab.