


# MICROBIOLOGY SWAB COLLECTION CHART

## BACTERIAL/FUNGAL COLLECTIONS



COPAN Dual Culture Swab System with Liquid Stuart's Medium

## CHLAMYDIA AND GONORRHOEAE COLLECTIONS BY AMPLIFIED DETECTION




APTIMA Unisex Swab Specimen Collection Kits for Endocervical and Male Urethral Swab Specimens



APTIMA Vaginal Swab Specimen Collection Kit (vaginal swabs)

## VIRAL COLLECTIONS



Rapid Flu – Puritan Sterile Foam-Tipped Applicator Swab



Rapid RSV – Puritan Sterile Foam-Tipped Applicator Swab



M4 Collection Kit Swab – Fisherbrand Polyester-Tipped Applicator Swab



M4 Collection Kit Swab – COPAN Nylon Flocked Mini-Tipped Flexible Plastic Applicator Swab

For more information regarding microbiology swab collections, please call (800) 934-6575 or visit [NationwideChildrens.org/Laboratory-Services](http://NationwideChildrens.org/Laboratory-Services).